# **MEDME**

# April 19, 2009

CGcount Determining the count of CpG dinucleotides for a set of genomic locations

## **Description**

The count of CpGs is determined in each window of size wsize, with or withouth weighting, for each probe according to its position, chromosome and genome realease

## Usage

```
CGcount(data, wsize = 1000, wFunction = "linear")
```

## **Arguments**

data An object of class MEDMEset

wsize number; the size of the smoothing window, in bp

wFunction string; the type of weighting function, to choose among linear, exp, log or none

## **Details**

Only human and mouse are currently supported. The respective genomic sequence metadata library needs to be downloaded from the Bioconductor website, installed and loaded (around 800Mb). Please note that only the last genome release should be used. LiftOver UCSC tool could be used for batch conversion of old genomic position to the last genome release.

#### Value

An object of class MEDMEset is returned where the count of CpGs for each probe has been saved on the CGcount slot.

## See Also

smooth

2 MEDME

#### **Examples**

```
data(testMEDMEset)
## just an example with the first 1000 probes
testMEDMEset = smooth(data = testMEDMEset[1:1000,])
library(BSgenome.Hsapiens.UCSC.hg18)
testMEDMEset = CGcount(data = testMEDMEset)
```

MEDME Determining the logistic model of MeDIP enrichment in respect to the

expected DNA methylation level

## **Description**

Probe-level MeDIP weighted enrichment is compared to the expected DNA methyllation level. The former is determined applying MeDIP protocol to a fully methylated DNA. The latter is determined as the count of CpGs for each probe. This is assumed to be the methylation level of each probe in a fully methylated sample.

## Usage

```
MEDME (data, sample, CGcountThr = 1, figName = NULL)
```

## **Arguments**

data An object of class MEDMEset

sample Integer; the number of the sample to be used to fit the model, based on the order

of samples in the smoothed slot

CGcountThr number; the threshold to avoid modelling probes with really low methylation

level, i.e. CpG count

figName string; the name of the file reporting the model fitting

## **Details**

The model should be applied on calibration data containing MeDIP enrichment of fully methylated DNA, most likely artificially generated (see references). Nevertheless, in case chromosome or genome-wide human tiling arrays are used a regular sample could be used too. In fact, human genomic DNA is known to be hyper-methylated but in the promoter regions. Of course the performance of the method is expected to be somehow affected by this approximation.

#### Value

The logistic model as returned from the multdrc function from the drc R library

## References

```
http://genome.cshlp.org/cgi/content/abstract/gr.080721.108v1
```

## See Also

```
smooth, CGcount
```

MEDME.readFiles 3

#### **Examples**

```
data(testMEDMEset)
## just an example with the first 1000 probes
testMEDMEset = smooth(data = testMEDMEset[1:1000, ])
library(BSgenome.Hsapiens.UCSC.hg18)
testMEDMEset = CGcount(data = testMEDMEset)
MEDMEmodel = MEDME(data = testMEDMEset, sample = 1, CGcountThr = 1, figName = NULL)
```

MEDME.readFiles

reading sgr or gff files for MEDME

## **Description**

allows to read sgr or gff files before submitting the data to MEDME analysis

## Usage

```
MEDME.readFiles(path = getwd(), files = NULL, format, organism)
```

## **Arguments**

path	string; the path where the files are stored; the current working directory is the default
files	vector; optional vector of file names
format	string; either sgr or gff to indicate the respective file formats
organism	string; either hsa or mmu for homo sapiens and mus musculus respectively

## **Details**

In case of GFF files (recommendend), tab-delimited files with header are expected with following fields: chromosome, probe ids, start and stop chromosomal positions, and score are expected in columns 1, 3, 4, 5 and 6 repectively. Multiple files are also expected to be in the same order of rows.

In case of sgr files (GFF is the preferred format), tab-delimited files with no header and chr, chr positions and score are expected in columns 1, 2 and 3 repectively. Multiple files are also expected to be in the same order of rows.

## Value

An object of class MEDMEset. The column headers in the logR slot are determined from the file names. in case of SGR files the are not probe names and progressive numbers are used in place of them. In case of GFF files the probe names are determined from the 3rd column.

4 MEDMEset-class

MEDME.writeFiles writeFiles sgr or gff files from MEDME output

## **Description**

allows to write sgr or gff files after MEDME analysis

## Usage

```
MEDME.writeFiles(data, output, path = getwd(), format, featureLength = NULL)
```

## **Arguments**

data An object of class MEDMEset

output string; the name of the data slot to be written on the disk, either logR, smoothed,

AMS or RMS

path string; the path where the files are stored; the current working directory is the

default

format string; either sgr or gff to indicate the respective file formats

featureLength

integer; in case of GFF file format the length of the features has to be provided

to determine start and end positions

## **Details**

One GFF or SGR file is provided for each sample of the data MEDMEset object.

In case of GFF files, tab-delimited files with header are provided with following fields for each probe: chromosome, empty field, probe ids, start and stop chromosomal positions, and score and empty fields.

In case of sgr files, tab-delimited files with no header and chr, chr positions and score are provided.

MEDMEset-class Class "MEDMEset"

## **Description**

This class is used in MEDME library to store MeDIP derived DNA-methylation estimates and to save further elaboration of these, in association with chromosomal and positional probe information

#### **Objects from the Class**

Objects can be created by calls of the form <code>new("MEDMEset", ...)</code>. This object could initially host the MeDIP normalized logRatio data, as returned by the <code>MEDME.readFiles</code> function. Afterwards, the same obejct is returned by most of the MEDME library function. Each time, a new slot is filled with additional data, as smoothed logR or Absolute/Relative Methylation Scores (AMS and RMS respectively). At the end of the analysis, usually after a call to the <code>predict</code> function, the <code>MEDME.writeFiles</code> function can be used to generate SGR or GFF files from this object.

MEDMEset-class 5

#### **Slots**

```
chr: Object of class "character": the probe-level chromosome asignments
pos: Object of class "numeric": the probe-level genomic position
logR: Object of class "matrix": the probe-level un-trasformed normalized MeDIP logRatios
    for each sample
smoothed: Object of class "matrix": the probe-level smoothed MeDIP logRatios for each
    sample

AMS: Object of class "matrix": the probe-level Absolute Methylation Score for each sample

RMS: Object of class "matrix": the probe-level Relative Methylation Score for each sample

CGcounts: Object of class "numeric": the probe-level count of CpGs

organism: Object of class "character": the organism that the probe genomic positions are
    referring to, either hsa or mmu for homo sapiens or mus musculus respectively
```

#### Methods

```
[ signature(x = "MEDMEset"): subsets the object based on its probes and/or samples
AMS signature (object = "MEDMEset"): extracts the Absolute Methylation Score from
    the AMS slot
CG signature (object = "MEDMEset"): extracts the probe CpG count from the CG-
    counts slot
chr signature(object = "MEDMEset"): extracts the probe chromosomal assignment
org signature(object = "MEDMEset"): extracts the organism
initialize signature(.Object = "MEDMEset"): automatically generates smoothed, AMS
    and RMS matrix when only the logR slot is filled
logR signature (object = "MEDMEset"): extracts the matrix of MeDIP un-transformed
    logRatios
pos signature(object = "MEDMEset"): extracts the probe genomic position
RMS signature (object = "MEDMEset"): extracts the Relative Methylation Score from
    the RMS slot
show signature (object = "MEDMEset"): prints a summary of the object content
smoothed signature(object = "MEDMEset"): extracts the Absolute Methylation Score
    from the AMS slot
```

## Author(s)

Mattia Pelizzola

#### References

```
http://genome.cshlp.org/cgi/content/abstract/gr.080721.108v1
```

## See Also

```
MEDME.readFiles, MEDME.writeFiles
```

## **Examples**

```
showClass("MEDMEset")
```

6 predict

predict

Applying the logistic model on MeDIP enrichment data

## **Description**

This allows the probe-level determination of MeDIP smoothed data, as well as absolute and relative methylation levels (AMS and RMS respectively)

## Usage

```
predict(data, MEDMEfit, MEDMEextremes = c(1,32), wsize = 1000, wFunction='linear
```

string; the type of weighting function, to choose among linear, exp, log or none

## **Arguments**

data An object of class MEDMEset

MEDMEfit the model obtained from the MEDME.model function

MEDMEextremes

vector; the background and saturation values as determined by the fitting of the model on the calibration data

wsize number; the size of the smoothing window, in bp

# Value

wFunction

An object of class MEDMEset. The resulting smoothed data, the absolute and relative methylation score (AMS and RMS) are saved in the smoothed, AMS and RMS slots, respectively.

## References

```
http://genome.cshlp.org/cgi/content/abstract/gr.080721.108v1
```

## See Also

```
smooth, CGcount, MEDME
```

## **Examples**

```
data(testMEDMEset)
## just an example with the first 1000 probes
testMEDMEset = smooth(data = testMEDMEset[1:1000, ])
library(BSgenome.Hsapiens.UCSC.hg18)
testMEDMEset = CGcount(data = testMEDMEset)
MEDMEmodel = MEDME(data = testMEDMEset, sample = 1, CGcountThr = 1, figName = NULL)
testMEDMEset = predict(data = testMEDMEset, MEDMEfit = MEDMEmodel, MEDMEextremes = c(1,32)
```

smooth 7

smooth	Determining weighted MeDIP data
--------	---------------------------------

## Description

MeDIP data from tiling arrays are smoothed by determining for each probe i the weighted average of the probes within a window of size wsize centered at i

## Usage

```
smooth(data, wsize=1000, wFunction='linear')
```

## **Arguments**

data An object of class MEDMEset

wsize number; the size of the smoothing window, in bp

wFunction string; the type of weighting function, to choose among linear, exp, log or none

#### Details

The un-smoothed data are read from the slot logR of the data MEDMEset and the resulting smoothed data are saved on the smoothed slot.

#### Value

An object of class MEDMEset. In particular, the smoothed data are saved on the smoothed slot.

## **Examples**

```
data(testMEDMEset)
# just an example with the first 1000 probes
testMEDMEset = smooth(data = testMEDMEset[1:1000,])
```

testMEDMEset

Dataset of class MEDMEset for testing MEDME

## **Description**

This dataset contains a subset of the data reported in references. It contains normalized un-smoothed probe-level MeDIP enrichment for almost 50000 probes. This is a random subset of a custom Nimblegen chromosome X tiling array. It is a two channels array with an resolution of 100bp and oligos of 60nt. The M value is reported only. The fullyMet column of the logR slot contains data from a calibration experiments where MeDIP has been applied to a fully methylated sample. The last two columns NBMEL and YUSAC2 contain DNA-methylation experimental data for two cell strains: NBMEL are newborn normal melanocytes cells and YUSAC2 a melanoma strain. Data was processed with within and between array normalization. The full dataset contains almost 380K probes. See references for details. Chromosome, genomic position and logR of probes can be accessed with the methods chr, pos and logR respectively.

Please note that the original genomic coordinates were mapped to the hg17 human genome. These have been converted to hg18 using the LiftOver USCS tool available online for batch conversion.

8 testMEDMEset

# Usage

data(testMEDMEset)

# Format

**MEDMEset** 

## References

http://genome.cshlp.org/cgi/content/abstract/gr.080721.108v1

# **Index**

```
*Topic classes
                                           smooth, 1, 2, 6, 7
   MEDMEset-class, 4
                                           smoothed (MEDMEset-class), 4
[, MEDMEset, ANY, ANY, ANY-method
                                           smoothed, MEDMEset-method
       (MEDMEset-class), 4
                                                  (MEDMEset-class), 4
AMS (MEDMEset-class), 4
                                           testMEDMEset, 7
AMS, MEDMEset-method
       (MEDMEset-class), 4
CG (MEDMEset-class), 4
{\tt CG}, {\tt MEDMEset-method}
       (MEDMEset-class), 4
CGcount, 1, 2, 6
CGcounts (MEDMEset-class), 4
chr (MEDMEset-class), 4
chr, MEDMEset-method
       (MEDMEset-class), 4
initialize, MEDMEset-method
       (MEDMEset-class), 4
logR (MEDMEset-class), 4
logR, MEDMEset-method
       (MEDMEset-class), 4
MEDME, 2, 6
MEDME.readFiles, 3, 4, 5
MEDME.writeFiles, 4, 4, 5
MEDMEset-class, 4
org (MEDMEset-class), 4
org, MEDMEset-method
       (MEDMEset-class), 4
organism (MEDMEset-class), 4
pos (MEDMEset-class), 4
pos, MEDMEset-method
       (MEDMEset-class), 4
predict, 4, 6
RMS (MEDMEset-class), 4
RMS, MEDMEset-method
       (MEDMEset-class), 4
show, MEDMEset-method
```

(MEDMEset-class), 4